

# Hemato – Flow

## Novi pristup u analizi DKS-a

Anica Remenar  
Beckman Coulter

WBC Est		
Neut	52	60
Lymph	35	29
Mono	7	8
Eos	1	3
Baso	1	
NRBC		
Band		
Imm Gran		
Var Lymph		
Blast		

*Myelo. 4*

Abnormalni  
uzorci

Brojanje 100  
stanica pod  
mikroskopom

Pojavom  
sumnjivih  
stanica razmaz  
se šalje  
citologu

100 CELLS COUNT AT  
MICROSCOPE:

1<sup>ST</sup> RUN 3 EO, 0 MYELOCYTES

2<sup>ND</sup> RUN 1 EO, 4 MYELOCYTES

WBC	3.21	L	10 <sup>3</sup> /μL
NE %	19.88	L	%
LY %	57.77	H	%
MO %	12.87	H	%
EO %	8.53	H	%
BA %	0.95		%
NRBC %	0.00		%

# Ograničenja ručne analize diferencijalne kvne slike

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*Fig. 1: Rümke Table – shows the statistical inaccuracy which results when counting out a defined number (n). Thus, the number of e.g. band cells can vary in a slide between 4.9 and 17.6 when 100 cells are counted out and 10 band cells are actually present. (Variation of the number due to subjective assessment is not taken into account here.)*

# Automatizacija : Hemato-Flow

Analyser LH 750  
(or DxH800)

or

FC 500


REMISOL


Cytodiff

Autre

Sample Prep  
FP 1000

# Hemato-Flow Cyto Diff panel

1 bočica

6 monoklonskih antitijela

5 boja

16 populacija

- **CD16PC5** : zreli i nezreli granulociti, citotoksični limfociti, limfociti T/NK
- **CD2+CD294+PE** : eozinofili, bazofili, aktivirani T, T i NK stanice
- **CD45PC7** : leukociti
- **CD36FITC** : monociti
- **CD19ECD** : limfociti B i B prekursori.



Abnormalni  
uzorci

Analiza 20,000  
stanica

Analiza stanica  
protočnim  
citometrom

A: LEUKOCYTES IN COLORS WBC  
DIFFERENTIAL IN 30.000 CELLS


WBC	3.21	L	10 <sup>3</sup> /μL
NE %	19.88	L	%
LY %	57.77	H	%
MO %	12.87	H	%
EO %	8.53	H	%
BA %	0.95		%
NRBC %	0.00		%

CELL TYPE	FCM TYPE
Lymphocytes	
■	B-lymphocytes
■	Non cytotoxic T-lymphocytes
■	Cytotoxic T/NK lymphocytes
Granulocytes	
■	Neutrophils
■	Eosinophils
■	Basophils
Monocytes	
■	CD16neg monocytes
■	CD16pos monocytes
Immature cells	
■	Immature granulocytes
■	Myeloid precursors
■	B-cell precursors
■	T-cell precursors
■ + ■ + ■	Total circulating precursors

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


## Refining the White Blood Cell Differential: The First Flow Cytometry Routine Application

Mikael Roussel,<sup>1\*</sup> Cyrille Benard,<sup>1</sup> Beatrice Ly-Sunnaram,<sup>1</sup> Thierry Fest<sup>1,2\*</sup>

- 800 hematoloških uzoraka (1 dan) - normalno 25 % treba ručno diferenciranje (160 – 200 slideova dnevno)
- Nakon CytoDiff – samo 22 uzorka trebala ručni diff (+ 16 koji su imali RBC flags)

## “6 Markers/5 Colors” Extended White Blood Cell Differential by Flow Cytometry

Jean-Luc Faucher,<sup>1</sup> Charlotte Lacronique-Gazaille,<sup>1</sup> Elise Frébet,<sup>1</sup> Franck Trimoreau,<sup>1</sup> Magali Donnard,<sup>1</sup> Dominique Bordessoule,<sup>2</sup> Francis Lacombe,<sup>3</sup> Jean Feuillard<sup>1\*</sup>

### Iz zaključka:

- WBC differential by FCM is at least as reliable as by standard cytology
- FCM was better than standard cytology in detection and quantification of circulating blast cells or immature granulocytes, with a first lineage orientation in the former case

# Zaključak

- Istraživanja<sup>1,2</sup> pokazuju da je kvantitativna detekcija abnormalnosti diferencijalne krvne slike neutrofila, eozinofila, bazofila, monocita i limfocita metodom protočne citometrije jednako pouzdana kao i standardna mikroskopija krvnih razmaza
- Cyto Diff je bolji od citologije u detekciji nezrelih granulocita i blasta

- 1. Jean-Luc Faucher et al., „6 Markers/5 Colors” Extended White Blood Cell Differential by Flow Cytometry, Cytometry Part A, 71A: 934\_944, 2007
- 2. Sven Björnsson, Total Nucleated Cell Differential for Blood and Bone Marrow Using a Single Tube in a Five-Color Flow Cytometer, Cytometry Part B (Clinical Cytometry) 74B:91–103 (2008)